





tralian Government Initiative



phn

An Australian G

phn NORTH WESTERN

ment Initiative

MEL BOURNE



An Australian Government Initiative

phn SOUTH EASTERN

MELBOURNE

Sterilisation

2016

Acknowledgment

This webinar has been developed by Eastern Melbourne PHN on behalf of the Victorian PHN Alliance, which is the collective platform for the six PHNs in Victoria.

Eastern Melbourne PHN does not take responsibility arising from the use of, or reliance on, this webinar by a third party. Any such use or reliance is the sole responsibility of that party. This webinar does not constitute medical advice. If you require medical advice, please consult an appropriate medical professional.

Objectives

- Define the principles of sterilisation and be able to describe the sterilisation process undertaken in the practice
- Differentiate between standard and transmission based precautions and their application in the sterilisation process
- Recognise the importance of validation and calibration of sterilisation procedures and process to ensure effectiveness of the sterilisation process

Content

- Risk assessment
- Processing reusable medical devices and equipment
- Staff training
- Equipment processing area
- Reprocessing equipment (sterilisation process)
- What monitoring and maintenance is required?
- What is involved in the calibration and validation process

RACGP Infection prevention and control standards 5th edition



www.racgp.org.au

Healthy Profession. Healthy Australia.

What do we mean by sterilisation?



Risk assessment and processing of reusable instruments and equipment



Must have a designated person responsible for sterilisation

able Consequences of any injury or ha ticant Moderate Majo e.g. first aid Extrem rate High Extrem V High Extrem V Moderate High V Moderate High

ne = immediate action

Undertake risk assessment-Spaulding classification

 Destination
 District Control (University)

 Marcine State
 Marcine State

 Marcine State
 Marcine State

<


Be able to describe every aspect of instrument and equipment reprocessing

Staff training






Treatment room





Reprocessing equipment



Daily cleaning procedures

• Dirty instrument container



Spray bottles



Correct dilution



Preliminary cleaning

Removing gross soil can be achieved by one or in combination of the following:

- Wipe off gross soil at the point of use
- Rinse all used instruments under gently running warm water in 'dirty' sink/container
- Immersing items in the labelled 'dirty instruments container'

Instrument cleaning

- Preparation for manual cleaning
- Use of standard precautions
- Washing and rinsing instruments
- Visually inspect



Cleaning


• Disposing of water

• Brush care

Cleaning the area



Sterile barrier systems (packaging)










Sterile barrier systems



Loading the steriliser





Unloading the steriliser

All sterilisers should produce automatically recorded evidence that the cycle parameters have been met for every cycle.

✓ Check bags for moisture and integrity

✓ Check chemical indicators

✓ Check for print out for pass/fail

Documenting the cycle

Record:

- Cycle date
- Load number
- Load contents
- Identity of person who prepared the load
- Results of cycle monitoring
- Class 1 chemical indicator change
- Results of other indicators used
- Condition of packaging
- Signature of person releasing or rejecting the load

Documentation





	RISK A	SSESSMENT M	ATRIX		
SEVERITY	Catastrophic (1)	Critical (2)	Marginal (3)	Negligible (4)	
Frequent (A)	High	High	Serious	Medium	
Probable (B)	High	High	Serious	Medium Low Low Low	
Occasional (C)	High	Serious	Medium		
Remote (D)	Serious	Medium	Medium		
improbable (E)	Medium	Medium	Medium		
Eliminated (F)	-	Eim	nated		

Sterilisation log

STERILISER LOG

✓ if all parameters are met

Date	Load Number & date on each pack. (Label with a nontoxic, solvent based felt tipped marking pen prior to sterilisation)	Description of items	Pre- cleaning, cleaning, drying & packaging completed by.	Class 1 chemical indicator changed.	Time and temperature printout or Class 4, 5 or 6 chemical indicator.	Unloaded & packs checked for complete dryness & integrity by:	Name of person releasing load.	Reason if part/whole load failed & action taken.	Description of maintenance and repairs.

Steriliser maintenance

- Ensures sterility of the equipment processed
- Assists in the longevity of the steriliser
- Should be carried out according to manufacturers instructions
- Annual Calibration and Validation



Validating the sterilisation process

Validation involves:

- Reviewing policies and procedures
- Performing all the procedures
- Checking the efficacy of the results
- Recording the results process
- Must be performed at installation and then annually

Single-use items

Must not be reprocessed

• Single patient use equipment may be reprocessed for use on the same patient

• Have enough single-use items available

Summary

- Risk assessment and the formulation of policies and procedures
- Staff training and competency
- Sterilisation process
- Calibration and validation of the steriliser



Thank you for watching

